



Built for larger, heavier instruments, this enhanced SterilGARD® solution delivers expanded workspace, high load capacity, and easy access.

Combining a trusted design with cost-effective support for OEM and third-party integrations.




KEY FEATURES


- Expanded work area with extra-deep (4" and 10") and extra-tall (5") dimensions to accommodate larger instruments
- Generous access and reach with up to 28" work depth, 34" work height, and a 26.5" access opening for easier instrument installation and handling
- High load capacity, supporting up to 300 lbs on the work surface
- Proven, familiar SterilGARD® design – users already know and trust it, eliminating the learning curve of specialized platforms
- Cost-effective solution, leveraging a modified SterilGARD® rather than a fully custom unit

Special Applications

Thermo KingFisher Presto
Automated Purification System



Sartorius AMBR 15
Bioreactor



Hamilton Microlab Prep
Liquid Handler



SterilGARD® Custom Sizes For Lab Automation

PRODUCTS DESIGNED	INSTRUMENT
SG404-23" VS max opening with raised diffuser	Agilent Bravo SRT liquid handler
SG404-4XD	Aspect Biosystems RX1 bioprinter
SG404-4XD- 23" VS opening	Sartorius Fill-it Automated Cryovial Filling System
SG404-4XD-Rear Wall pocket	Thermo KingFisher Flex
SG404-10XD	Miltenyi AutoMACS Pro, Thermo KingFisher Flex, Agilent MultiFlo FX, BioSpa & Citation system, Formulatrix F.A.S.T liquid handler
SG404-10XD-5XT	Sartorius AMBR 15 Bioreactor (small footprint), Desktop Health 3D Bioplotter, Beckman Coulter Cytoflex LX Flow Cytometer, Agilent Bravo liquid handler, Eppendorf epMotion liquid handler, Opentrons OT-2 liquid handler, Hamilton Microlab Prep liquid handler, Bausch+Strobel manual vial filler or capper
SG504-10XD	Any instrument that fits inside the SG404-10XD
SG504-10XT	Bausch+Strobel manual vial filler or capper, Hudson Robotics Scorpion liquid handler, Agilent Bravo NGS (Bravo on risers), Hamilton Microlab Nimbus
SG604-4XD-23" VS opening	Sartorius Fill-it Automated Cryovial Filling System, Sartorius AMBR15 Bioreactor
SG604-10XD	Miltenyi AutoMACS Pro, Cytex Aurora Cell Analyzer
SG604-5XT	Sartorius Fill-it Automated Cryovial Filling System (tight fit), Agilent Bravo
SG604-4XD-5XT	Sartorius Fill-it Automated Cryovial Filling System (best fit)
SG604-10XD-5XT	Sartorius AMBR 15 Bioreactor, Desktop Health 3D Bioplotter, Bausch+Strobel manual vial filler or capper, Beckman Coulter Cytoflex LX Flow Cytometer, Agilent Bravo liquid handler, Eppendorf epMotion liquid handler, Beckman Coulter Biomek 4000 liquid handler



OEM Partnership Ready

With compatibility for third-party instruments.

Baker Custom integrated SterilGARD® for 3P Innovation Cryovial filler; not for sale through Baker.



[in](#) [ig](#) [f](#) [v](#) bakerco.com

Scan to see more about our lab automation solutions

