

# 10 Ways

the AniGARD® e3 can make  
your lab more productive.

A photo e-book brought to you by:

THE BAKER COMPANY

[www.bakerco.com/anigard10ways](http://www.bakerco.com/anigard10ways)

# #1

**FEATURE:**  
3 sided work area

**BENEFIT:**  
Permits access  
by multiple  
operators and  
facilitates team-  
based projects

**IMPORTANT NOTE:**  
HEPA-filtered air  
throughout the work  
surface provides ISO  
Class 4 protection,  
reducing every  
operator's exposure  
to allergens.



# #2

**FEATURE:**

14" high access opening to a flat, spacious work surface.

**BENEFIT:**

Allows cages in a variety of sizes to be pulled in and out quickly.

AniGARD®  
ENERGY EFFICIENT ENGINEERED



# #3

**FEATURE:**  
Hinged  
viewscreen  
opens to 21 3/4"

**BENEFIT:**  
Greater visibility  
and access to  
work area -  
especially when  
moving large items  
in and out



# #4

**FEATURE:**  
Slanted  
viewscreen on  
all 3 sides

**BENEFIT:**  
Provides an  
ergonomic and  
direct line-of-sight  
panoramic view of  
the work space.



# #5

## FEATURE:

Low noise levels - 59 dBA when measured according to NSF method.

## BENEFIT:

Reduces effect of ambient noise on animals and is more comfortable for users.



# #6

**FEATURE:**


Well-lit workspace with optimal distribution - average of 100 foot-candles at work surface.

**BENEFIT:**

Significantly reduces eye strain.

**IMPORTANT NOTE:**

T5 energy efficient bulbs produce limited heat and are easy to replace.



# #7

**FEATURE:**

12" knee space opening

**BENEFIT:**

Provides seated access for small procedures.



# #8

## FEATURE:

Electric hydraulic lift comes standard

## BENEFIT:

Allows a work surface adjustment range to accommodate various sized operators comfortably.

## IMPORTANT NOTE:

Work surface height adjusts from 32 ½" to 44" with 5" casters.



# #9

**FEATURE:**

25' retractable power cord within the base of unit.

**BENEFIT:**

Provides greater range of use without unplugging and location avoids hits when moving around.

**IMPORTANT NOTE:**

Speaking of electricity - StediFLOW™ motor reduces energy and power requirements.



# #10

## FEATURE:

Wide tread stainless steel casters (5" standard, 6" optional) and optional pull bar.

## BENEFIT:

Easy maneuverability and movement throughout the facility.

## IMPORTANT NOTE:

Rubber bumpers protect the unit and the lab during minor collisions.



**Thank you for your interest!**

**Ready to take the next step?  
Request a quote at [www.bakerco.com](http://www.bakerco.com)**

A photo flip book brought to you by:

**THE BAKER COMPANY**

[www.bakerco.com/anigard10ways](http://www.bakerco.com/anigard10ways)