

IsoGARD®

Model: CEIG14

Electrical

Total Amps: 3
Amps: 15
Volts: 115
Watts: 345
BTU/HR: 1176

Classification

NSF Listed: n/a
NSF Class: III
NSF Type: n/a
Old NSF Class:
Old NSF Type: n/a

Weight

Net Weight: 950
Ship Weight: 1150

Containment

Particulates: product, personnel, environment
Vapors and Gases: product, personnel, environment (if exhausted, tre

Performance

Face Velocity: N/A
Recirculated Air %: 0
Exhausted Air %: 100

Exhaust

Contaminated Pressure Plenum: NEGATIVE
Exhaust Source: direct
Exhaust Destination: OUTSIDE
Exhaust Connection: HARD
Exhaust CFM:
Inches Of Water (Column):

Filter

HEPA Supply Filter(s): ABOVE WORK AREA
HEPA Exhaust Filter: YES

Features

StediVOLT: N/A
UniPressure: N/A
Stand: Standard
Warranty: 12 months