

AniGARD® VF

The AniGARD® VF Animal Transfer Station provides a spacious, HEPA-filtered work area for cage changing and animal handling. This single-sided, vertical-flow clean bench ensures ISO Class 5 (Class 100) air quality while enhancing productivity with a 14-inch access opening, ergonomic design, and mobile hydraulic stand.



FEATURES

- **ISO Class 5 (Class 100) Airflow** – HEPA-filtered vertical airflow protects animals and ensures particulate control.
- **Spacious Work Area** – 14" high access opening and flat work surface allow for easy cage handling.
- **Ergonomic Design** – Slanted 10° viewscreen and bright workspace enhance comfort and visibility.
- **Flexible & Mobile** – Adjustable stand with casters allows the unit to fit through standard doorways for easy transport between rooms.
- **Easy Maintenance** – Front-access HEPA filter and pre-filter system extend filter life and simplify service.
- **Quiet Operation** – Low noise levels improve user comfort and minimize stress on animals.
- **Adjustable Work Surface** – Electric/hydraulic lift system allows height adjustment from 22" to 40½" for user comfort and ease of transportation.

The AniGARD® VF Animal Transfer Station provides a spacious, HEPA-filtered workspace for safe and efficient cage changing. Its ergonomic design and mobility features enhance user comfort and lab productivity.



Scan to see the Animal Clean Benches Brochure

APPLICATION

The AniGARD® VF is ideal for animal research facilities, providing a controlled, HEPA-filtered environment for cage changing, animal handling, and contamination-sensitive procedures. Its spacious design and ergonomic features enhance workflow efficiency while ensuring both animal protection and particulate control.



TECHNICAL SPECIFICATION

MODEL NUMBER	AG4VF402	AG5VF502	AG6VF602
Size (nominal)	4'	5'	6'
Exterior Dimensions w/ Stationary Stand (w x d x h [min to max])	64" x 32" x [89(1/16)" - 97(9/16)"]	76" x 32" x [89(1/16)" - 97(9/16)"]	88" x 32" x [89(1/16)" - 97(9/16)"]
Interior Dimensions (w x d x h)	46" x 25(1/8)" x 31(13/16)"	58" x 25(1/8)" x 31(13/16)"	70" x 25 1/8" x 31(13/16)"
Shipping Weight			
Benchtop model	465 lbs	634 lbs	694 lbs
With channel stand	637 lbs	706 lbs	776 lbs
Access Opening (working position)	14"	14"	14"
Access Opening (loading/cleaning position)	30(3/8)"	30(3/8)"	30(3/8)"
Electrical Requirements	115 V, AC	115 V, AC	115 V, AC
Amps / Breaker	16 / 20	16 / 20	16 / 20

General Specifications

- Interior workspace made from 16-gauge stainless steel.
- ISO Class 5 (Class 100) air cleanliness in the workspace.
- 16-gauge cold-rolled steel exterior construction with white powder coated finish.
- HEPA filters with 99.99% minimum efficiency in capturing 0.3µm micron particles.
- Separate switches for blower and light; if blower is not on, light switch will not work.
- Washable pre-filter extends HEPA filter life.

RESOURCES

Scan the codes to view our helpful resources to maximise the use of the AniGARD®.



Visit the AniGARD VF website



Download the AniGARD e3 Sell Sheet